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```
import gc
import pandas
import pygame
import helpers as hlp
import constants as const
import global_config as glc
import instructions as instr
from pylink import *

def do_trial(trial, trials_design, sub_setting, part_info, part_results_file, cycle):
    letter_result = []

    message = "record_status_message 'Ensayo %d '" % (trial)
    getEYELINK().sendCommand(message)

    msg = "TRIALID %s" % str(trial)
    getEYELINK().sendMessage(msg)

    trial_discs = ''.join(
        ["/data/images/", "discs_", part_info['setting'], "_trial_", str(trial), ".png"])
    getEYELINK().sendMessage("!V TRIAL_VAR image %s" % trial_discs)

    if trial in const.DRIFT_CHECKINGS:
        hlp.drift_checking(glc.SCREEN)

    getEYELINK().setOfflineMode()
    msecDelay(50)
    error = getEYELINK().startRecording(1, 1, 1)
    if error:
        return error
    gc.disable()

    beginRealTimeMode(100)
    try:
        getEYELINK().waitForBlockStart(100, 1, 0)
    except RuntimeError:

        if getLastError()[0] == 0:
            hlp.end_trial()
            print("ERROR: No link samples received!")
            return TRIAL_ERROR
        else:
            raise

    eye_used = getEYELINK().eyeAvailable()
    if eye_used == const.RIGHT_EYE:
        getEYELINK().sendMessage("EYE_USED 1 RIGHT")
    elif eye_used == const.LEFT_EYE or eye_used == const.BINOCULAR:
        getEYELINK().sendMessage("EYE_USED 0 LEFT")
        eye_used = const.LEFT_EYE
    else:
        print("Error in getting the eye information!")
        return TRIAL_ERROR
    getEYELINK().flushKeybuttons(0)
    buttons = (0, 0)

    hlp.show_fixation_cross_timed(glc.SCREEN)

    startTime = currentTime()
    getEYELINK().sendMessage("SYNCTIME %d" % (currentTime() - startTime))

    hlp.show_discs_timed(glc.SCREEN, trial, trials_design, part_info)

    if sub_setting == 'bloqueo':
        letter_result = hlp.show_letter_timed(glc.SCREEN, trial, trials_design, part_info)
    if sub_setting == 'dirigido':
        hlp.show_arrow_timed(glc.SCREEN, trial, trials_design, part_info)
    if sub_setting == 'paralelo':
        hlp.show_asterisk_timed(glc.SCREEN)

    ticks = hlp.show_color_wheel(glc.SCREEN, trial, trials_design)
```

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```
running = True
while running:
    for event in pygame.event.get():
        if event.type == pygame.MOUSEBUTTONDOWN:
            running = False
```

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```
import random
import pygame
import math
import constants as const
import global_config as glc
from pygame.math import Vector2
from colormath.color_conversions import convert_color
from colormath.color_objects import HSLColor, AdobeRGBColor
from pylink import *

def drift_checking(screen):
    pygame.mouse.set_visible(0)
    while 1:
        if not getEYELINK().isConnected():
            return ABORT_EXPT
        try:
            error = getEYELINK().doDriftCorrect(screen.get_rect().w // 2, screen.get_rect().h // 2, 1, 1)
            if error != 27:
                break
        except:
            getEYELINK().doTrackerSetup()
        getEYELINK().doTrackerSetup()
        getEYELINK().setOfflineMode()
        msecDelay(50)

def end_trial():
    endRealTimeMode()
    pumpDelay(100)
    getEYELINK().stopRecording()
    while getEYELINK().getkey(): pass

def rgb2hsl(r, g, b):
    rgb = AdobeRGBColor(r, g, b)
    hsl_angle_color = convert_color(rgb, HSLColor).get_value_tuple()
    return hsl_angle_color[0]

def hsl2rgb(h, s=1, l=0.5):
    hsl = HSLColor(h, s, l)
    color = convert_color(hsl, AdobeRGBColor)
    rgb_color = AdobeRGBColor.get_upscaled_value_tuple(color)
    return rgb_color

def get_sub_setting(cycle, setting):
    if cycle in [1, 2]:
        sub_setting = const.setting_dic[setting][0]
    else:
        sub_setting = const.setting_dic[setting][1]
    return sub_setting

def rotate_color_wheel(screen, wheel):
    rotated_wheel = pygame.transform.rotate(wheel, random.randint(0, 360))
    rotated_wheel_rect = rotated_wheel.get_rect()
    rotated_wheel_rect.center = screen.get_rect().center
    return rotated_wheel, rotated_wheel_rect

def show_color_wheel(screen, trial, trials_design):
    position = Vector2()
    rotated_wheel, rotated_wheel_rect = rotate_color_wheel(screen, glc.WHEEL)
    screen.blit(rotated_wheel, (
        screen.get_rect().size[0] / 2 - rotated_wheel_rect.size[0] / 2,
        screen.get_rect().size[1] / 2 - rotated_wheel_rect.size[1] / 2))
    screen.blit(glc.QUESTION, glc.QUESTION_RECT)
    position.from_polar((const.DISTANCE, int(trials_design[trial][11][-3:])))
    pos = (int(screen.get_rect().center[0] + position[0]),
           int(screen.get_rect().center[1] - position[1]))
    pygame.gfxdraw.filled_circle(screen, pos[0], pos[1], const.RADIUS, const.BLACK)
    pygame.gfxdraw.aacircle(screen, pos[0], pos[1], const.RADIUS, const.WHITE)
    pygame.display.flip()
    getEYELINK().sendMessage("colwheel_start")
    col_wheel = ''.join(['../data/images/', 'col_wheel.png'])
    getEYELINK().sendMessage("IV IMGLOAD FILL %" % col_wheel)
    pygame.mouse.set_visible(1)
    ticks = pygame.time.get_ticks()
    return ticks

def show_color_wheel_fixed_colors(counter, event):
    instr_position = Vector2()
    pygame.mouse.set_visible(1)
    ins_wheel = pygame.transform.rotate(glc.WHEEL, 180)
    ins_wheel_rect = ins_wheel.get_rect()
    ins_wheel_rect.center = glc.SCREEN.get_rect().center
    glc.SCREEN.blit(ins_wheel, (
        glc.SCREEN.get_rect().size[0] / 2 - ins_wheel_rect.size[0] / 2, int(
            glc.SCREEN.get_rect().size[1] / 2 - ins_wheel_rect.size[1] / 2)))
    glc.SCREEN.blit(glc.QUESTION, glc.QUESTION_RECT)
    instr_position.from_polar((const.DISTANCE, 30))
    pos = (int(glc.SCREEN.get_rect().center[0] + instr_position[0]),
           int(glc.SCREEN.get_rect().center[1] + -instr_position[1]))
    pygame.gfxdraw.filled_circle(
```

```

    glc.SCREEN, pos[0], pos[1], const.RADIUS, const.BLACK)
pygame.gfxdraw.aacircle(
    glc.SCREEN, pos[0], pos[1], const.RADIUS, const.WHITE)
x, y = event.pos
dist = abs(math.sqrt(
    (x - glc.SCREEN.get_rect().center[0]) ** 2 +
    (y - glc.SCREEN.get_rect().center[1]) ** 2))
if glc.INNER_CIRCLE_RADIUS < dist < glc.OUTER_CIRCLE_RADIUS:
    pygame.mouse.set_visible(0)
    glc.SCREEN.fill(const.BLACK)
    counter += 1
return counter

```

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```

def wait_click(screen, trial, trials_design, part_info, part_results_file, letter_result, ticks):
    running = True
    while running:
        for event in pygame.event.get():
            if event.type == pygame.MOUSEBUTTONDOWN:
                x, y = pygame.mouse.get_pos()
                dist = abs(math.sqrt(
                    (x - screen.get_rect().center[0]) ** 2 +
                    (y - screen.get_rect().center[1]) ** 2))
                if glc.INNER_CIRCLE_RADIUS < dist < glc.OUTER_CIRCLE_RADIUS:
                    cw_response_time = pygame.time.get_ticks() - ticks
                    pygame.mouse.set_pos(screen.get_rect().center)
                    getEYELINK().sendMessage("colwheel_end")
                    pygame.mouse.set_visible(0)
                    hsl_color, target_color = compute_color(
                        screen, trial, trials_design, x, y)
                    transformed_error, non_transformed_error = compute_transformed_error(
                        hsl_color, target_color)
                    write_csv_part_results(
                        part_results_file, part_info, trial, trials_design,
                        hsl_color, target_color, non_transformed_error,
                        transformed_error, letter_result, cw_response_time)
                    return 0, transformed_error

```

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```

def show_discs_timed(screen, trial, trials_design, part_info):
    screen.fill(const.BLACK)
    show_discs(screen, trial, trials_design, part_info)
    pygame.time.delay(const.DISCSTIME)
    screen.fill(const.BLACK)
    getEYELINK().sendMessage("discs_end")

```

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```

def check_letter_response(trial, trials_design):
    correct_letter_type = trials_design[trial][8]
    ticks = pygame.time.get_ticks()
    letter_response = 0
    elapsed_time = 0
    pressed_click = 0
    while elapsed_time < const.LETTERTIME:
        elapsed_time = pygame.time.get_ticks() - ticks
        for event in pygame.event.get():
            if event.type == pygame.MOUSEBUTTONDOWN and event.button in [const.LEFT, const.RIGHT] and not pressed_click:
                letter_response = pygame.time.get_ticks() - ticks
                if event.button == const.LEFT:
                    letter_type = const.VOCAL
                if event.button == const.RIGHT:
                    letter_type = const.CONSONANT
                if letter_type == correct_letter_type:
                    letter_accuracy = 1
                else:
                    letter_accuracy = 0
            pressed_click = 1
    if letter_response == 0:
        letter_response = const.NO_RESPONSE_LETTER
        letter_accuracy = 0
        letter_type = const.NO_RESPONSE_LETTER
    return [letter_response, letter_type, letter_accuracy]

```

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```

def show_letter_feedback(trial, letter_response, letter_accuracy):
    if trial in const.FEEDBACK_TRIALS:

```

name (participant name)	<code>part_info['name'],</code>
par_id (3 digits ID, 001-999)	<code>part_info['par_id'],</code>
sex (masculino, femenino)	<code>part_info['sex'],</code>
age (18-99)	<code>part_info['age'],</code>
session (1)	<code>part_info['session'],</code>
setting (1-4)	<code>const.setting_id[part_info['setting']],</code>
hsl_color (0-360)	<code>hsl_color,</code>
target_color (0-360)	<code>target_color,</code>
non_transformed_error (0-360)	<code>non_transformed_error,</code>
transformed_error (0-180)	<code>transformed_error,</code>
colorwheel_response_time (milisegundos)	<code>str(cw_response_time),</code>
letter_presented (aA-zZ, excepto ñ)	<code>trials_design[trial][7],</code>
letter_type_presented (no_letter, consonant, vocal)	<code>trials_design[trial][8],</code>
letter_response_time (milisegundos)	<code>str(get_letter_result(trial, trials_design, letter_result)[0]),</code>
letter_type_selected (-999, consonant, vocal)	<code>str(get_letter_result(trial, trials_design, letter_result)[1]),</code>
letter_accuracy (0-1)	<code>str(get_letter_result(trial, trials_design, letter_result)[2]))</code>

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```
import sys
import pygame
import pygame.gfxdraw
import helpers as hlp
import constants as const
import global_config as glc
from pygame.math import Vector2

def present_instructions(cycle, setting, part_info):
    sub_setting = hlp.get_sub_setting(cycle, setting)
    present_pre_cycle_instructions(cycle)
    if cycle in [1,3]:
        present_sub_setting_instructions(sub_setting, part_info)

def import_instructions_file(INSTR_PATH):
    with open(INSTR_PATH, encoding = 'latin-1') as f:
        list_instr = f.read().splitlines()
    instructions = {}
    for b in list_instr:
        i = b.split(';')
        instructions[i[0]] = glc.FONT.render(i[1], True, const.WHITE)
    return instructions

def present_pre_cycle_instructions(cycle):
    instructions = import_instructions_file(const.PRE_CYCLE_INSTR_PATH[cycle])
    running = True
    counter = 0
    unders_instr = False
    pygame.mouse.set_visible(0)
    if cycle != 1:
        hlp.print_centered_msg(const.CLICK_MSG)

    while running:
        for event in pygame.event.get():
            if event.type == pygame.KEYDOWN and event.key == pygame.K_ESCAPE:
                running = False
                pygame.quit()
                sys.exit()
            if event.type == pygame.MOUSEBUTTONDOWN and not unders_instr and counter < 4:
                glc.SCREEN.fill(const.BLACK)
                if counter == 0:
                    hlp.print_msg_L1(instructions['1.0'])
                    hlp.print_msg_D(instructions['0.0'])
                if counter == 1:
                    hlp.print_msg_L1(instructions['2.0'])
                    hlp.print_msg_D(instructions['0.0'])
                if counter == 2:
                    hlp.print_msg_L1(instructions['3.0'])
                    hlp.print_msg_D(instructions['0.0'])
                if counter == 3:
                    unders_instr = True
                    running = False
                    counter += 1
                if cycle in [2,4] and counter == 2:
                    counter = 3

        pygame.display.flip()
    if cycle in [2,4]: msg_next_stage = const.CLICK_BEGIN_MSG
    else: msg_next_stage = const.CLICK_MSG
    next_stage = glc.FONT.render(msg_next_stage, True, const.WHITE)
    next_stage_rect = next_stage.get_rect()
    next_stage_rect.center = glc.SCREEN.get_rect().center
    glc.SCREEN.blit(next_stage, next_stage_rect)
    pygame.display.flip()
    if cycle in [2,4]:
        running = True
        while running:
            for event in pygame.event.get():
                if event.type == pygame.KEYDOWN and event.key == pygame.K_ESCAPE:
                    running = False
                    pygame.quit()
                    sys.exit()
                if event.type == pygame.MOUSEBUTTONDOWN:
                    running = False

def present_sub_setting_instructions(sub_setting, part_info):
    instructions = import_instructions_file(const.INSTR_PATH[sub_setting])
    running = True
    counter = 0
    unders_instr = False
    instr_presented = False
    pygame.mouse.set_visible(0)
    while running:
        for event in pygame.event.get():
            if event.type == pygame.KEYDOWN and event.key == pygame.K_ESCAPE:
                running = False
                pygame.quit()
```

```

sys.exit()
if event.type == pygame.MOUSEBUTTONDOWN and not unders_instr and counter < 8:
    glc.SCREEN.fill(const.BLACK)
    if counter == 0:
        hlp.print_msg_L1(instructions['0.1'])
        hlp.print_msg_D(instructions['0.0'])
    if counter == 1:
        hlp.show_fixation_cross(glc.SCREEN)
        hlp.print_msg_L1(instructions['1.0'])
        hlp.print_msg_D(instructions['0.0'])
    if counter == 2:
        hlp.show_discs_fixed_colors(glc.SCREEN)
        hlp.print_msg_L1(instructions['2.0'])
        hlp.print_msg_D(instructions['0.0'])
    if counter == 3:
        hlp.print_msg_L1(instructions['3.0'])
        hlp.print_msg_L2(instructions['3.1'])
        hlp.print_msg_D(instructions['0.0'])
        if sub_setting == 'bloqueo':
            hlp.show_letter(glc.SCREEN, 0, 'A', 'left', part_info)
            hlp.print_msg_L3(instructions['3.2'])
        if sub_setting == 'paralelo':
            hlp.show_asterisk(glc.SCREEN)
        if sub_setting == 'dirigido':
            hlp.show_arrow(glc.SCREEN, 0, '330', part_info)
            hlp.print_msg_L3(instructions['3.2'])
            hlp.print_msg_L4(instructions['3.3'])
    if counter == 4:
        hlp.print_msg_L1(instructions['4.0'])
        hlp.print_msg_D(instructions['0.0'])
    if counter == 5:
        hlp.print_msg_L1(instructions['5.0'])
        hlp.print_msg_D(instructions['0.0'])
    if counter == 6:
        hlp.show_color_wheel_fixed_colors(counter, event)
    if counter == 7:
        hlp.print_msg_L1(instructions['6.0'])
        instr_presented = True
        counter += 1
if (event.type == pygame.KEYDOWN and event.key == pygame.K_n and
    instr_presented):
    pygame.mouse.set_visible(0)
    glc.SCREEN.fill(const.BLACK)
    hlp.print_msg_L1(instructions['7.0'])
    pygame.display.flip()
    counter = 0
    unders_instr = False
    instr_presented = False
if (event.type == pygame.KEYDOWN and event.key == pygame.K_s and instr_presented):
    unders_instr = True
    running = False
pygame.display.flip()

```

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